BARD® BIOPTY-CUT®

Disposable Core Biopsy Needle

Instructions for Use

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BARD
BIOPSY SYSTEMS

Instructions for Use:

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

A. General Information and Device Description:

The BARD® BIOPTY-CUT® Biopsy Needle is a single patient use core biopsy needle designed exclusively for use with the BARD® BIOPTY® Instrument. The needle may have centimeter markings for depth measurement and is available in several needle gauge sizes and lengths. The BARD® BIOPTY-CUT® Biopsy Needle can be purchased with an attached spacer. The spacer is useful in maintaining the distance between the two needle hubs promoting ease of handling and simplifying targeting.

B. How Supplied

1. The needle is supplied sterile and non-pyrogenic unless the package has been opened or damaged. Sterilized using Ethylene Oxide. **For single use only. Do not reuse. Do not resterilize.**
2. The needle is available in various gauge sizes and lengths.

C. Indications for Use:

The core needle biopsy device is intended for use in obtaining biopsies from soft tissues such as liver, kidney, prostate, spleen, lymph nodes and various soft tissue tumors. It is not intended for use in bone.

D. Contraindications:

Good medical judgment should be exercised in considering biopsy on patients who are receiving anticoagulant therapy or who have a bleeding disorder.

E. Warnings:

1. **After placing the BARD® BIOPTY-CUT® Biopsy needle into the BARD® BIOPTY® Instrument, care should be taken to prevent premature activation. Leaving the safety pin in the "Ready" position could result in premature activation of the device and possible injury to the patient and/or user.**
2. **The BARD® BIOPTY-CUT® Biopsy Needle has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.**
3. **Do not resterilize the BARD® BIOPTY-CUT® Biopsy Needle. After resterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.**

NOTE: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and with applicable local, state, and federal laws and regulations.

F. Precautions:

1. Use only BARD® BIOPTY-CUT® Biopsy Needles with the BARD® BIOPTY® Instrument. We cannot recommend the use of biopsy needles manufactured by other manufacturers.
2. This product should be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of core needle biopsy, in particular, those relating to the specific organ being biopsied.
3. The introduction of the needle into the body should be carried out under imaging control (Ultrasound, X-Ray, CT, etc.).
4. Never test the instrument while the needle is installed in the instrument. This could result in needle damage or injury to the patient and/or user.
5. Unusual force applied to the stylet or unusual resistance against the stylet while extended out of the supportive cannula may cause the stylet to bend at the specimen notch. A bent specimen notch may interfere with needle function.

G. Potential Complications:

Potential complications of core biopsy are site specific and may consist of hematoma; hemorrhage; infection; adjacent tissue injury; pain; bleeding; hemoptysis; hemothorax; non-target tissue, organ or vessel perforation; and air embolism. Air embolism is a rare but serious potential complication of lung biopsy procedures. Rapid deterioration of neurological status and/or cardiac arrhythmia may be indicative of air embolism. Prompt diagnosis and treatment must be considered if the patient exhibits signs or symptoms of air embolism.

H. Equipment Required:

- Appropriate imaging modality accessories
- Surgical gloves and drapes
- Local anesthetic as needed
- BARD® TRUGUIDE® Coaxial cannula (optional)
- Scalpel
- Sample collection container
- Other equipment as necessary

I. Directions for Use:

BARD® BIOPTY-CUT® Biopsy Needle Preparation:

1. Prior to use, determine the appropriate gauge and length of the needle required for the specific biopsy to be performed. Using aseptic technique, remove the needle from the package and from the protective sheath.
2. **Precaution:** Before using, inspect the needle for damaged point, bent shaft or other imperfections that would prevent proper function. If the needle is damaged or bent, **DO NOT USE**.
3. Before loading needle into the instrument, slide cannula over the stylet several times noting the proper free movement. **DO NOT USE** if needle specimen notch does not pass into the outer cannula freely.

NOTE: A slight drag may be felt if the inner stylet is removed completely from the outer cannula. This is not a defect. The tip is made with a slight back bend to make sure the notch clears the outer cannula.

Biopsy Procedure:

1. Skin Preparation:

- A. The biopsy procedure must be performed using appropriate aseptic technique.
- B. Prepare site as required. Adequate anesthesia should be administered prior to incision of the skin.
- C. **RECOMMENDATION:** For ease of insertion, puncture the skin with a scalpel at the entry site.

2. Positioning:

A. With BARD® BIOPTY® Instrument and BARD® BIOPTY-CUT® Biopsy Needle (spacer optional):

Open cover of energized (cocked) instrument. Install needle by placing the outer cannula hub in the forward needle hub carriage and the inner stylet hub in the rear needle hub carriage. See Figure B. Then remove the spacer (if used) by grasping the tab and pulling straight out. See Figure C. Close cover. Introduce needle, under imaging control where appropriate, through incision until the needle point is proximal to the area to be biopsied.

-OR-

B. With Only the BARD® BIOPTY-CUT® Biopsy Needle with Spacer:

While maintaining the stylet and cannula hub orientation with the spacer, introduce needle, under imaging control where appropriate, through incision until the needle point is proximal to the area to be biopsied. Open cover of energized (cocked) instrument. Install needle by placing the outer cannula hub in the forward needle hub carriage and the inner stylet hub in the rear needle hub carriage, being careful to maintain the needle's position and orientation. See Figure B. Remove the spacer by grasping the tab and pulling straight out. See Figure C. Close cover.

FIGURE A


FIGURE B



FIGURE C

3. Activation:

While maintaining the instrument's position and needle's orientation, move safety pin toward "ready" position. Depress the actuator button to cause both the stylet and cannula to automatically advance.

4. Completion:

- A. Withdraw the instrument and needle from the patient. Open the instrument and remove the needle.
- B. While holding the stylet hub stationary, slide the cannula back to expose the biopsy section.
- C. The collection of multiple needle cores may help to ensure the detection of any cancer tissue. A "negative" biopsy in the presence of suspicious radiographic findings does not preclude the presence of carcinoma.

NOTE: If collecting multiple samples, inspect the needle for damaged point, bent shaft or other imperfections after each sample is collected. **DO NOT USE** needle if any imperfection is noted.

NOTE: Post-biopsy patient care may vary with the biopsy technique utilized and the individual patient's physiological condition. Observation of vital signs and other precautions should be taken to avoid and/or treat potential complications that may be associated with biopsy procedures.

J. Cleaning and Lubrication of the BARD® BIOPTY® Instrument:

Precautions:

- **DO NOT** use 2% glutaraldehyde or other bactericidal solutions, as it may cause the instrument not to function properly.
- **DO NOT** sterilize by irradiation.
- **DO NOT** clean with concentrated chlorhexidine gluconate or any other agent containing alcohol.

Cleaning:

Inspect the instrument for signs of deterioration or damage (cracking, blistering, separation of coating, pitting). **DO NOT USE** if deterioration or damage is observed. Thoroughly wash the instrument in a mild, warm detergent. Then rinse all parts well with hot water. Thoroughly dry the product inside and out.

Lubrication:

Use quality medical grade instrument lubrication. Bard recommends that the instrument be cleaned and lubricated after every use.

Warranty:

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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Copyright ©2012 C. R. Bard, Inc. All rights reserved.
Copyright ©2012 C. R. Bard, Inc. Tous droits réservés.
Copyright ©2012 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright ©2012 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright ©2012, C. R. Bard, Inc. Reservados todos los derechos.
Copyright ©2012 C. R. Bard, Inc. Alle rechten voorbehouden.